

Stereo Indirect System

High Resolution , Wide-Angle. Surgical observation field without deformation ,
To make the posterior surgery perfect.



Technical Parameters

Lens	Field of view	Image magnification	Lens diameter
Wide-angle lens	F175 : 102° F200 : 120°	0.43X	19mm
Ultra small diameter wide-angle lens	F175 : 95° F200 : 112°	0.42X	13mm
Mid view lens	F175 : 80° F200 : 95°	0.74X	19mm

Compatible microscope brand and model

ZEISS/TOPCON/EDER SOI4	LEICA SOI5	Möller SOI6	Others SOI7
---------------------------	---------------	----------------	----------------

Wide-angle lens



Ultra small diameter wide-angle lens



Mid view lens



Forte Biyoteknoloji San. Ve Tic. Ltd. Şti.

www.fortebiyoteknoloji.com

Adalet Mah. Manas Bulv. Folkart Towers A Kule 47/B K:26 D:2601 PK:35530 Bayraklı İzmir / Türkiye
Tel: +90.232.332.04.82
Fax: +90.232.332.04.90
Email: info@fortebiyoteknoloji.com



Note

For your better use of the instrument, please be sure to read the instructions carefully before use.